



# eResept 2.3.101 (#4686)

## Release Notes


27.4.2012

---

Copyright ©2012 Pula ehf. All rights reserved.

## No component defined

---

 It is now possible to configure the printer configuration in two ways:

- 1) Per computer, suitable where a computer is connected to certain printers.
- 2) Per computer and user, suitable where users are connecting via remote desktop solutions. Each user then has their own setting based on what computer they are logged onto.


The default value is "per computer" so installations that will not require user to have their own printer configuration do not need to do any changes to the configuration.

The configuration is visible under the advanced view in the system configuration part in the admin client.

[ RENO-3813 ]

 All types of inbox notifications can now be flagged as handled directly in the doctor's inbox.

[ RENO-3820 ]

 When a doctor signs and sends, all the inbox notifications related to draft prescriptions created by an assistant are automatically marked as read.

[ RENO-3800 ]